

HOW TO MEASURE TRIPODS

TUBE LENGTH



WARNING

Length of the tube at constant diameter
(do not measure tightening ring or tube base)

TUBE PERIMETER

- 1/ cut a strip of paper of 15 cm long and 1 cm wide, and roll it tight around the tube
- 2/ make a mark at the junction of each paper strip tips (fig1)
- 3/ cut the paper strip at the mark
- 4/verify each tip of the paper strip are edge to edge when rolled around the tube and are not overlapping (fig 2)
- 5/ measure the new length of the paper strip which gives you the perimeter of the tube

Figure 1




Figure 2

